

## Nitrosamine Analysis in Tobacco and Tobacco Product using LC-MS/MS

Okta Lian Atikah, S.Si | Field Application Support LC-MS/MS



18 February 2020.



## Outline

- Nitrosamine and Tobacco-specific Nitrosamine
- Determination of tobacco-specific nitrosamines in tobacco and tobacco products by LC-MS/MS (CRM No. 72)
- Other contaminants analysis in Tobacco using LC-MS/MS







## Nitrosamine

- **Nitrosamines** are organic compounds of the chemical structure R2N–N=O
- Most nitrosamines are carcinogenic.
- Nitrosamines can be found in water, foods, cosmetics, tobacco, packing materials and drugs.
- Tobacco-specific nitrosamines (TSNAs) are nitrosamine compound group that only found in tobacco products, and possibly in some other nicotine-containing products



Nitrosamine Structure





## Tobacco-Specific Nitrosamine (TSNAs)

- Tobacco-specific nitrosamines (TSNAs) are formed from nicotine and related compounds by a nitrosation reaction that occurs during the curing and processing of tobacco
- The IARC (International Agency for Research on Cancer) has identified 8 tobacco-specific nitrosamines in tobacco and tobacco smoke.
- Two of them have been classified as Group 1 carcinogens:
  - N-nitrosonornicotine (NNN)
  - 4-methyl-N-nitrosamino-1-(3-pyridyl)-1-butanone (NNK)



 In January 2017, US FDA (Food and Drug Administration) has proposed rule to establish limit of N-Nitrosonornicotine (NNN) in finished smokeless tobacco product

How to determine TSNAs in Tobacco and Tobacco Product?



nswers for Science nowledge for Life



## CORESTA Recommended Method

CORESTA (Cooperation Centre for Scientific Research Relative to Tobacco) has released 3 recommended methods (CRM) for determination of TSNAs

- CORESTA Recommended Method (CRM) No. 63: Determination Of Tobacco Specific Nitrosamines in Cigarette Mainstream Smoke – GC-TEA Method (Newest version released January 2019)
- CORESTA Recommended Method (CRM) No. 75: Determination Of Tobacco Specific Nitrosamines in Cigarette Mainstream Smoke by LC-MS/MS (Newest version released August 2019)
  - The reproducibility data was better for LC-MS/MS than for GC-TEA
- CORESTA Recommended Method (CRM) No. 72: Determination Of Tobacco Specific Nitrosamines in Tobacco and Tobacco Products by LC-MS/MS (Newest version released July 2017)

\*coresta.org





# CORESTA Recommended Method No.72

Determination of tobacco-specific nitrosamines (TSNAs) in tobacco and tobacco products by LC-MS/MS

- The TSNAs determined with this method are
  - N-nitrosonornicotine (NNN)
  - 4-(N-methylnitrosamino)-1-(3-pyridyl)-1-butanone (NNK)
  - N-nitrosoanatabine (NAT)
  - N-nitrosoanabasine (NAB).





Answers for Science. Knowledge for Life.™



## Method Workflow – Targeted

Determination of tobacco-specific nitrosamines (TSNAs) in tobacco and tobacco products by LC-MS/MS



Sample and Standard Preparation

Sample Analysis Using LC-MS/MS

Data Process





## Calibration Standard Preparation

- Stock solution for standard and IS is diluted in ACN ٠
- Mixed TSNA std. Nr.2: The concentration is approximately 400 ng/ml of NNN, ٠ NNK, NAT and 100 ng/ml of NAB in ACN:Water (30:70)
- Internal std spiking solution 2000 ng/mL diluted in ACN. ٠

∣Cal. Std.	Volume of Mixed TSNA std. Nr. 2	Volume of Internal std spiking solution 2000 ng/ml	Volume acetonitrile	lume Conc. onitrile NNN		Conc. NAT	Conc. NAB	
	(ml)	(ml)	(ml)	(ng/ml)	(ng/ml)	(ng/ml)	(ng/ml)	
Cal 1	0,125	1,00	22	0,5	0,5	0,5	0,125	
Cal 2	0,250	1,00	22	1,0	1,0	1,0	0,250	
Cal 3	0,50	1,00	22	2,0	2,0	2,0	0,50	
Cal 4	1,00	1,00	22	4,0	4,0	4,0	1,00	
Cal 5	2,00	2,00 1,00		8,0	8,0	8,0	2,00	
Cal 6	5,00	5,00 1,00		20	20	20	5,00	
Cal 7	25,0 1,00		15	100	100	100	25,0	

Each of the cal std is diluted to volume 100 mL with 100 mM ammonium ٠





## Sample Preparation

Samples	Preparation
Tobacco and tobacco products	Tobacco and tobacco products are grounded Grinding procedure shall not generate heat or cause sample degradation
Smokeless tobacco products	Shall be analyzed together with their pouch (paper) and shall be cut into two halves directly into the extraction flask
Cigar and Cigarette Filler	Cigar and cigarette filler are grounded. Testing may also involve the analysis of the entire cigar where the wrapper and filler are ground together.

\*The test samples shall be stored protected from light





## Sample Extraction



## LC Conditions

- Mobile Phase:
  - Mobile phase A: Water
  - Mobile phase B: 0,1 % (v/v) acetic acid in methanol
- Column: C18 HPLC Column, 2,5 μm particle size, 2,1 mm × 50 mm, or equivalent
- Column Temperature: 60,0 °C
- Injection Volume: 10  $\mu L$
- LC Gradient

Time (Min)	Flow rate (mL/min)	Conc. A (%)	Conc. B (%)
0	0.22	100	0
3.0	0.22	10	90
4.0	0.22	10	90
5.0	0.22	0	100
6.0	0.22	100	0
10.0	0.22	100	0





## MS conditions – MRM Scan Mode

			· · · ·
Compound Name	Quantitation Transition (m/z)	Qualification Transition (m/z)	Internal Standard Reference
NNK	208 > 122	208 > 79	NNK-d4
NNK-d4	212 > 126	n/a	n/a
NNN	178 > 148	178 > 105	NNN-d4
NNN-d4	182 > 152	n/a	n/a
NAT	190 > 160	190 > 79	NAT-d4
NAT-d4	194 > 164	n/a	n/a
NAB	192 > 162	192 > 133	NAB-d4
NAB-d4	196 > 166	n/a	n/a

- Quantitation Transition will be used for Quantitation analysis in other to determine concentration of TSNAs in samples
- Qualification Transition will be used for confirmation using ion ratio





# US FDA – Determination NNN in Smokeless Tobacco and Tobacco Filler

# US FDA has adapted CRM No.72 with some modifications to determine NNN in smokeless tobacco and tobacco filler by LC-MS/MS

The same extraction process based on CRM No. 72 with modifications to the calibration ranges and instrumental approach.

LIB# 4623 Page 1 of 9

FDA/ORA/ORS

LABORATORY INFORMATION BULLETIN

Determination of N-nitrosonornicotine (NNN) in Smokeless Tobacco and Tobacco Filler by  ${\rm HPLC-MS/MS^1}$ 

John A. Mathis, Ph.D., Rayman A. Stanelle, Ph.D., and Jean-Marie D. Dimandja, Ph.D. U.S. Food and Drug Administration, Office of Regulatory Affairs, Southeast Regional Laboratory, Tobacco Branch, 60 8<sup>th</sup> St. Atlanta, GA 30309. John.Mathis@fda.hhs.gov

Note: The Laboratory Information Bulletin is a tool for the rapid dissemination of laboratory methods (or information) which appears to work. It does not report complete scientific work. The user must assure himself/herself by appropriate validation procedures that LIB methods and techniques are reliable and accurate for his/her intended use. Reference to any commercial materials, trade names, equipment, or procedure does not necessarily constitute approval, endorsement, or recommendation by the Food and Drug Administration or the U.S. Department of Health and Human Services.



### SCIEX 6500 QTRAP



Answers for Science. Knowledge for Life." PT. LABORINDO SARANA

## Confirmation using ion ratio



LABORINDO SARANA

SCIEX Answers for Science. Knowledge for Life."

## SCIEX OS – Data Processing

[MQ4] Untitled Method	X	Build Processing
Workflow	For each component, configure the parameters to optimize peak integration	Method
Components	NNN 1 Integration NNN 1 (261.2 / 158.1) from Sample A (Carisprodol validasi.wiff (sample 23)) Area: 3,740e5. Heinht 9,477e4. RT: 3,85 min	Components
Library Search	Apply peak parameters to all of the components Minimum Peak Width 3 points 8a4	<ul> <li>Integration</li> </ul>
Calculated Columns	Minimum Peak Height 100.00 7e4 -	Ion Ratio
Flagging Rules	Gaussian Smooth Width 0.0 points ≩ 5e4 -	acceptance
Advanced	Noise Percentage 40.0 % <u>E</u> Con Baseline Subtract Window 2.00 min 3e4 -	Accuracy
Non-targeted Peaks	Peak Splitting     2     points       Retention Time (RT)     0e0	acceptance
	1 2 3 4 5 6 7 8 9 10 11 12 Apply Time, min	• etc

### **Result Table**

- Area
- Retention time
- Calculated Concentration
- Ion Ratio
- Ion Ratio Acceptance
- Calibration Curve



Answers for Science. Knowledge for Life.™

ን	28	rows Filters: 0	Qualify for R	ules 🔀 🕻		°c "C	ılı. C	C,,H, 🗾 👻	071:		More		• 👹	$\left  \times \right $
1	ndex	Sample Name 🛛	Component Name ⊽	Component Type ⊽	Actual Concentr ⊽	Expected RT ▽	Area ⊽	, Retent ⊽ Time	Calculated Concentrat ▽	Accuracy N	7 Ion Ratio ⊽	lon Ratio	Co Ac ⊽	^
	22	QC 0.2ppb	NNN 2	Qualifiers	0.20	3.85	6.568e5	3.86	2.055e-1	02.73	0.8497	~		
	23	QC 0.2ppb	NNN 1	Quantifiers	0.20	3.85	7.776e5	3.85	2.084e-1	04.22	0.8504	<ul> <li>Image: A set of the set of the</li></ul>		
	24	QC 0.2ppb	NNN 2	Qualifiers	0.20	3.85	6.613e5	3.85	2.070e-1	103.51	0.8504	<ul> <li>Image: A second s</li></ul>		
►	25	Sample A	NNN 1	Quantifiers	N/A	3.85	3.740e5	3.85	8.707e-2	N/A	0.8377	$\sim$		
	26	Sample A	NNN 2	Qualifiers	N/A	3.85	3.133e5	3.85	8.356e-2	N/A	0.8377	<ul> <li>Image: A set of the set of the</li></ul>		
	27	Sample B	NNN 1	Quantifiers	N/A	3.85	2.214e6	3.86	6.405e-1	N/A	0.8398	<ul> <li>Image: A second s</li></ul>		
	28	Sample B	NNN 2	Qualifiers	N/A	3.85	1.860e6	3.86	6.323e-1	N/A	0.8398	<ul> <li>Image: A second s</li></ul>		
									•					~



# Which other contaminant can be analyzed using LC-MS/MS?

### Nicotine analysis in tobacco using LC-MS/MS

STUDIA UNIVERSITATIS BABEŞ-BOLYAI, PHYSICA, L, 4b, 2005

#### DETERMINATION OF NICOTINE FROM TOBACCO BY LC-MS-MS

Laurian Vlase<sup>1</sup>, Lorena Filip<sup>2</sup>, Ioana Mîndruțău<sup>3</sup>, Sorin E. Leucuța<sup>1</sup>

Faculty of Pharmacy, University of Medicine and Pharmacy 'I uliu Hatieganu'', Cluj-Napoca <sup>1</sup>Department of Pharmaceutical Technology <sup>2</sup>Department of Environmental Chemistry <sup>3</sup>Department of Physical Chemistry

#### ABSTRACT

A new high performance liquid chromatography coupled with mass spectrometry method (LC/MS/MS) for quantification of nicotine from tobacco was elaborated. It was utilized an Atlantis HILIC, 100 mm x 3.0 mm i.d., 3  $\mu$ m column with a mobile phase containing acetonitrile/ solution 0.2% formic acid in water. The ionization was optimized using ESI(+) and enhanced selectivity was achieved using tandem mass spectrometric analysis. The precursor to product ion transitions of m/z 163 -> (105.8, 131.8) were used to measure the nicotine concentrations. The quantification was made using the external standard method. The calibration curves were made on range 0.04-4  $\mu$ g/ml. Nicotine content was determined in 40 brands of cigarettes available in Romania. With our methodology, we obtained values of nicotine in tobacco between 7.5 and



Chromatogram result of Nicotine in Tobacco Leaf Using LC-MS/MS

ABORINDO SARANA



Answers for Science. Knowledge for Life.™

# Which other contaminant can be analyzed using LC-MS/MS?

### Pesticide Analysis using LC-MS/MS

# SCIEX LC-MS/MS allow us to analyze hundred of analytes in one run



#### Article

#### Analysis of 118 Pesticides in Tobacco after Extraction With the Modified QuEChRS Method by LC-MS-MS

Fei Yang<sup>1</sup>, Zhaoyang Bian<sup>1,2\*</sup>, Xiaoshui Chen<sup>1</sup>, Sansan Liu<sup>1</sup>, Yang Liu<sup>1</sup> and Gangling Tang<sup>1</sup>

<sup>1</sup>China National Tobacco Quality Supervision & Test Center, Zhengzhou 450000, PR China, and <sup>2</sup>No.2, Fengyang Street, Zhengzhou New & Hi-Tech Industries Development Zone, Henan, China

\*Author to whom correspondence should be addressed. Email: zybian@126.com

Received 3 February 2013; revised 10 June 2013

A liquid chromatography-tandem quadrupole mass spectrometry (LC-MS-MS) multi-residue method for the simultaneous target analysis of a wide range of pesticides in tobacco has been developed. Gradient elution has been used in conjunction with positive mode electrospray ionization tandem mass spectrometry to detect up to 118 pesticides in tobacco. The recoveries obtained for each pesticide ranged between 70 and 118% at two spiked concentration levels. Good linear relationships were observed with correlation coefficients  $r^2 > 0.992$  for all analytes. The established method was successfully applied to the determination of pesticide residues in real tobacco samples in order to validate the suitability for routine analysis. based on a liquid partitioning with acctonitrile followed by a dispersive solid-phase extraction (d-SPE) clean-up with primary secondary amine (P8A). Modifications to the original method to ensure efficient extraction of pH-dependent compounds (by using different buffer solutions) (14, 15) or addition of water to dry samples in order to obtain the necessary moisture (16, 17) have been introduced. The QuEChERS method is particularly popular for determination of polar, middle polar and non-polar pesticide residues in various food matrices (17–21) because of its simplicity, inexpensiveness, amenability to high throughput and relatively high-efficiency results with a minimal number of steps. In this work a simpla liquid chromosprenethy-trandem under











# Scheduled MRM<sup>™</sup> Algorithm



- Adjusts detection windows automatically depending on retention time
- Optimizes dwell times for each analyte and cycle time
  - Allows detecting many more MRM transitions
  - Allows using faster LC
  - Gives best S/N, accuracy and reproducibility





## Related method reference

#### Nitrosamines Analysis

- Mathis, et al. Determination of N-nitrosonornicotine (NNN) in Smokeless Tobacco and Tobacco Filler by HPLCMS/MS. (US FDA)
- Li, et al. 2015. Simultaneous Determination of Alkaloids and Their Related Tobacco-Specific Nitrosamines in Tobacco Leaves Using LC–MS-MS. Yunnan Academy of Tobacco Agricultural Sciences
- Guo, et al. 2018. Rapid Analysis of Genotoxic Nitrosamines by HPLC-MS/MS. (SCIEX)

#### Pesticide Analysis

- CORESTA Guide N° 5 Technical Guide for Pesticide Residues Analysis on Tobacco and Tobacco Products
- Development of a Fast and Cost-Effective Multi-Residue Method to Determine Pesticides in Tobacco by LC/MS/MS
- Dr. Hans Rainer Wollseifen. Determination of Pesticides in Tobacco Automated with FREESTYLE QuECHERS and LC-MS/MS.

#### **Nicotine Analysis**

 McGuffrey, et al. 2014. Validation of a LC-MS/MS Method for Quantifying Urinary Nicotine, Six Nicotine Metabolites and the Minor Tobacco Alkaloids—Anatabine and Anabasine—in Smokers' Urine





## Conclusion

- Determination of TSNAs can be done using LC-MS/MS according to CRM No. 72 and has adapted by US FDA
- Data processing can be done using SCIEX OS Software
- Other potent carcinogenic contaminants, such as nicotine and pesticide can also analyzed by LC-MS/MS
- ScheduleMRM allow us to analyze hundred of analyte in one run





## Reference

- CORESTA Recommended Method No. 72 Version 4. July 2017
- Mathis, et al. Determination of N-nitrosonornicotine (NNN) in Smokeless Tobacco and Tobacco Filler by HPLCMS/MS. US FDA
- Vlase, et al. 2005. Determination of Nicotine from Tobacco Using LC-MS/MS
- Yang, et al. 2013. Analysis of 118 Pesticides in Tobacco after Extraction With the Modified QuEChRS Method by LC-MS/MS



Answers for Science. Knowledge for Life.<sup>™</sup>

## THANK YOU FOR LISTENING!

It's Time to See the Future Differently

Email: okta.atikah@laborindo.com

Learn more about LC-MS/MS through www.sciex.com